




Ultraviolet-Visible Spectroscopy






□ The Spectroscopic Process

1. In **UV-vis spectroscopy**, the sample is **irradiated** with the broad spectrum of the **UV-vis radiation**
2. If a particular **electronic transition** matches the **energy** of a certain **band**, it will be **absorbed**
3. The remaining **light** passes through the sample
4. From this residual radiation **a spectrum** is obtained with “**gaps**” at these **discrete energies** – this is called **an absorption spectrum**.



□ Observed electronic transitions

From the **molecular orbital diagram**, there are several possible **electronic transitions** that can occur, each of a **different relative energy**



Transitions

$\sigma \rightarrow \sigma^*$

UV photon required, high energy

Methane at 125 nm (CH_4)

Ethane at 135 nm (C_2H_6)

$n \rightarrow \sigma^*$

Saturated compounds with unshared e^-

Absorption between 150 nm to 250 nm

ϵ between 100 and 3000 $\text{L cm}^{-1} \text{mol}^{-1}$

Shifts to shorter wavelengths with polar solvents

Minimum accessibility

Halogens, N, O, S

$n \rightarrow \pi^*$, $\pi \rightarrow \pi^*$

Organic compounds, wavelengths 200 to 700 nm

Requires unsaturated groups

$n \rightarrow \pi^*$ low ϵ (10 to 100)

$\pi \rightarrow \pi^*$ higher ϵ (1000 to 10000)

INSTRUMENTATION



Spectrometer: An instrument used for measuring transmittance or absorbance of a sample as function of the wavelength of the electromagnetic radiation

Components of a spectrophotometer

Source: It generates electromagnetic radiation


Dispersion device: It select the wavelength required from the broad band of radiation source

Sample area: Where the sample is kept

Detector: One or more detectors that measure the intensity of the radiation


SPECTROMETER DESIGN

Conventional Spectrometer




- Polychromatic source of light falls on an **entrance slit** which transmits a **narrow band** of light.
- The light then pass through the **sample** to a detector.
- The detector measures the **absorbance** of the sample **by comparing the light that reaches the detector from the sample and the blank (only solvent)**

Diode array Spectrometer



- Polychromatic source of light falls on a the sample, the transmitted radiation pass through an **entrance slit** of the dispersion device.
- The detector measures the **absorbance** of the sample **by comparing the light that reaches the detector from the sample and the blank (only solvent)**

DUAL BEAM SPECTROMETER



Radiation Sources


Two sources are required to scan the entire UV-Visible band:

- **Deuterium arc** lamp– covers the **UV – 200-330** (yields a good intensity continuum in the UV region)
- **Tungsten-halogen** lamp – covers **330-700** (yields good intensity over part of the UV spectrum and over the entire visible range)
- An alternate light source: **Xenon** lamp
 - Pros:** Yields a good continuum over the entire **UV** and **visible** regions.
 - Cons:** **High noise** from currently available **Xenon lamps** compared to **deuterium** or **tungsten** lamps

□ **Monochromator:** consists of an **entrance** slit, a **dispersion** device, and an **exit** slit.



Dispersion devices

- ❖ Dispersion devices cause different wavelengths of light to be dispersed at different angles.
- ❖ When combined with an appropriate exit slit, these devices can be used to select a particular wavelength of light from a continuous source.
- ❖ Two types of commonly used dispersion devices:
 1. Prisms
 2. Holographic gratings



Sample array

1. Sample cells can be made of **plastic, glass or quartz**
2. **Glass** absorbs strongly below **320 nm**
3. The cells **lowest in cost** are made of **plastic**, usually an **acrylic**. These cells **are not resistant** to all solvents and **absorb strongly** below **300 nm**
4. Only **quartz** is transparent in the full **200-700 nm** range; **plastic** and **glass** are only suitable for **visible regions spectra**



DETECTORS


❖ A **detector** converts a **light signal** into an **electrical signal**. It gives a **linear response** over a wide range **with low noise** and **high sensitivity**.

❖ Spectrophotometers normally contain


Photomultiplier tube detector

or a

Photodiode detector.



Photomultiplier tube



Photodiode


Principle of UV-Visible spectroscopy


- **Beer-Lambert's law**

$$A = \epsilon c l$$

- A = Absorbance of the sample = $\log_{10}(I_0/I)$
- I = Intensity of the radiation emerging from the sample
- I_0 = Intensity of the radiation incident on the sample
- ϵ = extinction coefficient or molar absorptivity in $M^{-1}cm^{-1}$
- c = concentration of the sample in moles/litr
- l = length of the light path through the sample in cm

UV/Vis of Acetone






Applications of UV-Vis spectroscopy

Terminology


- **Chromophore:** The **group of atoms** producing an **absorption** is called a 'Chromophore'. Eg. **C=O**, **C=C** etc.
- **Auxochrome:** Groups that helps in **extending the conjugation** by means of the **lone pairs** present. Eg. **NH₂**, **OH**, **OR**, **Cl**, **Br** etc.
- **Hypsochromic shift:** When the λ_{max} of an absorption shift to **shorter** wavelength
- **Bathochromic shift:** When the λ_{max} of an absorption shift to **longer** wavelength.
- **Hyperchromic shift:** When the **absorbance intensity** of λ_{max} shift to **higher** values
- **Hypochromic shift:** When the **absorbance intensity** of λ_{max} shift to **lower** values



Descriptive term	Nature of the shift
Bathochromic shift (Red shift)	Towards longer wavelength
Hypsochromic shift (Blue shift)	Towards shorter wavelength
Hyperchromic effect	Towards higher absorbance
Hypochromic effect	Towards lower absorbance



Applications of UV-Vis spectroscopy

Effect of conjugation



Applications of UV-Vis spectroscopy

Effect of conjugation





- ✓ Wavelength of absorption shift to longer wavelength with increase in conjugation

Applications of UV-Vis spectroscopy

Effect of substituents



AUXOCHROME

e.g. Benzene $\lambda_{\max} = 255 \text{ nm}$





Toluene $\lambda_{\max} = 261 \text{ nm}$

Phenol $\lambda_{\max} = 270 \text{ nm}$





Phenoxide ion $\lambda_{\max} = 287 \text{ nm}$

Aniline $\lambda_{\max} = 280 \text{ nm}$



Anilinium ion $\lambda_{\max} = 254 \text{ nm}$

Chlorobenzene $\lambda_{\max} = 265 \text{ nm}$



Naphthalene $\lambda_{\max} = 312 \text{ nm}$


- ❖ The **substituents** with **lone pair of electrons** may undergo **conjugation** with **π -electrons** of the phenyl ring causing '*Bathochromic shift*'
- ❖ **Blocking** of the **non-bonding pair of electrons** by **protonation** cause '*Hypsochromic shift*'

Applications of UV-Visible spectroscopy


Conjugation in 'enones'

Conjugation of C=O with C-C π bonds results in absorbance at higher wavelengths

Individually, C-C pi bonds and C-O pi bonds each group absorb at < 200 nm ($\pi \rightarrow \pi^*$)




absorption $\lambda_{\max} < 200$ nm



absorption $\lambda_{\max} < 200$ nm (C=O $\pi \rightarrow \pi^*$)


270 nm (C=O $n \rightarrow \pi^*$)
(weak)

When combined in conjugation, absorbance moves to longer wavelengths
(228 nm for mesityl oxide, below)



Mesityl oxide
absorption λ_{\max} 228 nm
($\pi \rightarrow \pi^*$)

Similar to butadiene



absorption λ_{\max} 217 nm

❖ Exercises

- Q1: A 2.5×10^{-4} M solution of a substance in a 1 cm length cell at λ_{max} 245 nm has absorbance 1.17. Calculate molar extinction coefficient for this transition.
- Q2: A 0.01 M solution of a compound transmits 20% of the radiation in a container with a path length equal to 1.5 cm. Calculate molar extinction coefficient of the compound.
- Q3: How do you differentiate between Dual beam, conventional, diode array Spectrophotometers?
- Q4: What do you understand by Chromophore and auxochrome? Illustrate with one example each.
- Q5: What do you understand by bathochromic and hypsochromic shift? Illustrate with one example each.
- Q6: Account for the high molar extinction coefficient observed in β -carotene, lycopene and retinal?
- Q7: What are different electronic transition possible in Benzophenone, KMnO_4 ?

